

RP00966

Leader in Biomolecular Solutions for Life Science



# Recombinant Human Thrombospondin-1/THBS1 Protein

Catalog No.: RP00966

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
Human	7057	P07996

### Tags

C-10×His

### Synonyms

THBS1;THBS;THBS-1;TSP;TSP-1;TSP1

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

< 1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.2 μm filtered solution of PBS, pH7.4. Contact us for customized product form or formulation.

### Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

## Background

## Basic Information

### Description

Recombinant Human Thrombospondin-1/THBS1 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Asn19-Pro1170) of human Thrombospondin-1/THBS1 (Accession #P07996) fused with a 10×His tag at the C-terminus.

### Bio-Activity

### Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.


## Contact



[www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Recombinant Human  
Thrombospondin-1/THBS1 Protein was  
determined by SDS-PAGE under reducing  
conditions with Coomassie Blue.