

RP00977

Leader in Biomolecular Solutions for Life Science



# Recombinant Human IFN-lambda 1/IL-29 Protein

Catalog No.: RP00977

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
Human	282618	Q8IU54

### Tags

C-His

### Synonyms

IFNL1;IL-29;IL29

## Background

This protein is also known as cytokine Zcyto21, Interferon lambda-1, IFN-lambda-1, and IFNL1, is a secreted protein which belongs to the IL-28 / IL-29 family. IL-29 is a cytokine with immunomodulatory activity. IL-29 is highly similar in amino acid sequence to the IL-28. IL-28 and IL-29 are induced by viral infection and showed antiviral activity. IL-28 and IL-29 interacted with a heterodimeric class II cytokine receptor that consisted of IL-1 receptor beta (IL-1Rbeta) and an orphan class II receptor chain, designated IL-28Ralpha. IL-29 plays an important role in host defenses against microbes and its gene is highly upregulated in cells infected with viruses. IL-29 may play a role in antiviral immunity. IL-29 up-regulates MHC class I antigen expression. It is a Ligand for the heterodimeric class II cytokine receptor composed of IL1RB and IL28RA. The ligand / receptor complex seems to signal through the Jak-STAT pathway.

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

Please contact us for more information.

### Formulation

Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

### Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

## Basic Information

### Description

Recombinant Human IFN-lambda 1/IL-29 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Gly20-Thr200) of human IL-29 (Accession #NP\_742152.1) fused with a 6xHis tag at the C-terminus.

### Bio-Activity

### Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.


## Contact



[www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Recombinant Human IFN-lambda 1/IL-29  
Protein was determined by SDS-PAGE with  
Coomassie Blue, showing a band at 27-37  
kDa.