

RP01404

Leader in Biomolecular Solutions for Life Science



Recombinant Human FGF-3 Protein

Catalog No.: RP01404

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	2248	P11487

Tags

No tag

Synonyms

FGF3;HBGF-3;INT2

Background

Fibroblast Growth Factor 3 (FGF-3) belongs to the large FGF family which has at least 23 members. All FGF family members are heparin-binding growth factors with a core 120 amino acid (aa) FGF domain that allows for a common tertiary structure. FGFs are expressed during embryonic development and in restricted adult tissues. They act on cells of mesodermal and neuroectodermal origin to regulate diverse physiologic functions including angiogenesis, cell growth, pattern formation, embryonic development, metabolic regulation, cell migration, neurotrophic effects and tissue repair. Signaling receptors for FGFs are type I transmembrane receptor tyrosine kinases belonging to the Ig superfamily. Four distinct but related classes of FGF receptors, FGF R1, 2, 3, and 4, exist. Through alternative splicing, multiple isoforms for FGF R1, 2 and 3, with distinct ligand recognition profiles, are also generated.

Product Information

Source	Purification
<I>E. coli</I>	> 95% by SDS-PAGE.

Endotoxin

Please contact us for more information.

Formulation

Lyophilized from a 0.22 µm filtered solution of 25mM Tris, 1M NaCl, 1mM TCEP, pH 8.0

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

Basic Information

Description

Recombinant Human FGF-3 Protein is produced by <I>E. coli</I> expression system. The target protein is expressed with sequence (Asp28-Arg212[N-Met]) of human FGF-3 (Accession #NP_005238.1) fused with no additional amino acid.

Bio-Activity

Storage


Store the lyophilized protein at -20°C to -80°C for long term. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human FGF-3 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 24kDa.