

RP01441

Leader in Biomolecular Solutions for Life Science



# Recombinant Human Apolipoprotein A-II/APOA2 Protein

Catalog No.: RP01441

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
Human	336	P02652

### Tags

C-hFc&Avi

### Synonyms

APOA2; Apo-AII; ApoA-II; apoAII; ApoA-II; apoAII

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

<0.1EU/μg

### Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

### Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

## Background

Apolipoprotein A-II (Apo-AII for short), also known as Apolipoprotein A2, is a secreted protein which belongs to the apolipoprotein A2 family. It exists as a disulfide-linked homodimer; and also can form a disulfide-linked heterodimer with APOD. APOA2 is the 2nd most abundant protein of the high density lipoprotein particles. This protein may stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism. Defects in APOA2 gene might cause apolipoprotein A-II deficiency or hypercholesterolemia.

## Basic Information

### Description

Recombinant Human Apolipoprotein A-II/APOA2 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Ala19-Gln100) of human Apolipoprotein A-II/APOA2 (Accession #NP\_001634.1) fused with a Fc, Avi tag at the C-terminus.

### Bio-Activity

### Storage

Store the lyophilized protein at -20°C to -80°C for long term. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.


## Contact



[www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Recombinant Human Apolipoprotein A-II/APOA2 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 45-50kDa.