

RP01771

Leader in Biomolecular Solutions for Life Science



Recombinant Mouse Periostin/PN/POSTN Protein

Catalog No.: RP01771

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Mouse	50706	Q62009-2

Tags

C-His

Synonyms

PN; PLF; Osf2; OSF-2;
A630052E07Rik[POSTN]

Product Information

Source	Purification
HEK293 cells	

Endotoxin

< 0.001EU/μg

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the vial before opening.
Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.
Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

Background

Periostin (POSTN), also known as OSF2 (osteoblast specific factor 2), is a heterofunctional secreted extracellular matrix (ECM) protein comprised of four fasciclin domains that promotes cellular adhesion and movement, as well as collagen fibrillogenesis. Postn is frequently overexpressed in various types of human cancers, stimulating metastatic growth by promoting cancer cell survival, invasion and angiogenesis, and can be a useful marker to predict the behavior of cancer.

Basic Information

Description

Recombinant Mouse Periostin/PN/POSTN Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Asn24-Ile653) of mouse Periostin/PN/POSTN (Accession #NP_056599.1) fused with a 6xHis tag at the C-terminus.

Bio-Activity

Storage


Store the lyophilized protein at -20°C to -80°C for 12 months. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.
Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Mouse Periostin/PN/POSTN Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 75 kDa.