

# Recombinant *Schistosoma japonicum* GST Protein

Catalog No.: RP02000LQ **Recombinant**

## Sequence Information

Species	Gene ID	Swiss Prot
<i>Schistosoma japonicum</i>		P08515

### Tags

No tag

### Synonyms

GST; GST 26; Sj26 antigen; SjGST

## Product Information

Source	Purification
<i>E. coli</i>	> 95 % as determined by SDS-PAGE.

### Endotoxin

Please contact us for more information.

### Formulation

Supplied as 0.22um filtered solution in PBS, pH 7.4. Contact us for customized product form or formulation.

### Reconstitution

## Background

GST isoenzymes appear to play a central role in the parasite detoxification system. Other functions are also suspected including a role in increasing the solubility of haematin in the parasite gut.

## Basic Information

### Description

### Bio-Activity

### Storage

This product is stable at  $\leq -70^{\circ}\text{C}$  for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.


## Contact



| [www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Recombinant *Schistosoma japonicum* GST Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 29 kDa.