

RP02364

Leader in Biomolecular Solutions for Life Science



Biotinylated Recombinant Cynomolgus FcRn Protein

Catalog No.: RP02364

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Cynomolgus		Q8SPV9&Q8S PW0

Tags

C-His&Avi

Synonyms

FcRn alpha chain; alpha-chain; Fc fragment of IgG; receptor; transporter; alpha; FCGRT; FCRN; FCRNimmunoglobulin receptor; intestinal; heavy chain; IgG Fc fragment receptor transporter alpha chain; IgG receptor FcRn large subunit p51; major histocompatibil

Product Information

Source	Purification
HEK293 cells	> 95% by Tris-Bis PAGE; > 95% by SEC-HPLC

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Background

Basic Information

Description

Biotinylated Recombinant Cynomolgus FcRn Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Met1-Ser297 (FCGRT)&Met1-Met119 (B2M)) of Cynomolgus FcRn fused with a His tag and Avi tag at the C-terminal.

Bio-Activity

Storage


Store the lyophilized protein at -20°C to -80°C for long term. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

Contact




www.abclonal.com

Validation Data



Biotinylated Recombinant FcRn on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of Biotinylated FcRn is greater than 95% as determined by SEC-HPLC.