

RP02461

Leader in Biomolecular Solutions for Life Science



Biotinylated Recombinant Human Notch 1 Protein

Catalog No.: RP02461

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	4851	P46531

Tags

C-His&Avi

Synonyms

AOS5;AOVD1;TAN1;hN1;NOTCH1;Activated NOTCH1;notch 1; AOS5; AOVD1; TAN1; hN1; notch 1

Product Information

Source	Purification
HEK293 cells	> 95% by Tris-Bis PAGE;> 95% by SEC-HPLC

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Background

Basic Information

Description

Biotinylated Recombinant Human Notch 1 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ala19-Gln526) of Human Notch 1 fused with His tag and Avi tag at the C-terminal.

Bio-Activity

Storage


Store the lyophilized protein at -20°C to -80°C for long term. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

Contact




www.abclonal.com

Validation Data



Biotinylated Human Notch 1 on Tris-Bis PAGE under reduced conditions. The purity is greater than 95%.



The purity of Biotinylated Human Notch 1 is greater than 95% as determined by SEC-HPLC.