

RP03179

Leader in Biomolecular Solutions for Life Science



Recombinant Human Vimentin Protein

Catalog No.: RP03179

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	7431	P08670

Tags

C-His

Synonyms

CTRCT30; HEL113;Vimentin

Product Information

Source	Purification
Baculovirus-Insect Cells	> 95% by SDS-PAGE.

Endotoxin

< 1 EU/μg

Formulation

Lyophilized from sterile 40% acetonitrile, 0.1% TFA

Reconstitution

Centrifuge the vial before opening.
Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.
Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

Background

Vimentin is a type III intermediate filament (IF) protein found in various non-epithelial cells, especially mesenchymal cells. A vimentin monomer, has a central α -helical domain and carboxyl (tail) domains. Two monomers compose the basic subunit of vimentin assembly. Vimentin is crucial for supporting and anchoring the position of the organelles in the cytosol. Vimentin provided cells with a resilience absent from the microtubule or actin filament networks, when under mechanical stress in vivo. Therefore, in general, it is accepted that vimentin is the cytoskeletal component responsible for maintaining cell integrity. Vimentin is also responsible for stabilizing cytoskeletal interactions. It is found that vimentin control the transport of low-density lipoprotein. It has been used as a sarcoma tumor marker to identify mesenchyme.

Basic Information

Description

Recombinant Human Vimentin Protein is produced by Baculovirus-Insect Cells expression system. The target protein is expressed with sequence (Met1-Glu466) of human Vimentin (Accession #NP_003371) fused with a His tag at the C-terminus.

Bio-Activity

Storage


Store the lyophilized protein at -20°C to -80°C for long term. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human Vimentin Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 55-70 kDa.