

RP10000

Leader in Biomolecular Solutions for Life Science



Recombinant human 15N-labeled Ubiquitin protein

Catalog No.: RP10000

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human		

Tags

No Tag

Synonyms

Product Information

Source	Purification
--------	--------------

<I>E. coli</I>

Endotoxin

Formulation

Supplied in 20 mM Tris, 150 mM NaCl, 2 mM BME, 10% Glycerol. Contact us for customized product form or formulation.

Reconstitution

Background

15N-labeled ubiquitin can be used in NMR experiments for determining interactions between ubiquitin and other proteins. It can also be used as a standard for quantification of cellular ubiquitin levels in quantitative mass spectrometric assays.

Basic Information

Description

Bio-Activity

Storage

This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt.
For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature.


Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant human 15N-labeled Ubiquitin Protein was determined by SDS-PAGE with Coomassie Blue.