

RP10135

Leader in Biomolecular Solutions for Life Science



Recombinant human Heat shock 70 kDa protein 1A/HSP70 protein

Catalog No.: RP10135

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	3303	

Tags

No Tag

Synonyms

Background

HSP70 a critical chaperone protein that has a high affinity for unfolded polypeptide chains. It binds extended peptide segments with a net hydrophobic character exposed by polypeptides during translation and membrane translocation, or following stress-induced damage. In cooperation with other chaperones, Hsp70 stabilizes preexistent proteins against aggregation and mediates the folding of newly translated polypeptides in the cytosol as well as within organelles. Hsp70 and Hsp40 recruit some misfolded proteins to the ubiquitin E3 ligase CHIP for ubiquitination.

Basic Information

Description

Bio-Activity

Storage

This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.

Product Information

Source	Purification
<i>E. coli</i>	

Endotoxin

Formulation

Supplied in 20 mM Tris, 150 mM NaCl, 2 mM BME, 10% Glycerol. Contact us for customized product form or formulation.

Reconstitution

Contact



www.abclonal.com