

RP10219K

Leader in Biomolecular Solutions for Life Science



# Ubiquitin Linkage Screening Kit III

Catalog No.: RP10219K

## Sequence Information

Species	Gene ID	Swiss Prot
Human		
proteins		

### Tags

### Synonyms

## Product Information

Source	Purification
E. coli recombinant	95% based on SDS-PAGE

### Endotoxin

### Formulation

20 mM Tris, 150 mM NaCl, 2 mM  $\beta$ ME,  
10% Glycerol

### Reconstitution

## Background

Ubiquitin Linkage Screening Kit III is designed for two purposes: 1) to identify the ubiquitin lysine residue that is used to link the polyubiquitin chain conjugated on a protein substrate; 2) to determine whether a protein substrate is polyubiquitinated or conjugated with multiple monoubiquitin.

## Basic Information

### Description

### Bio-Activity

### Storage

This product is stable at  $\leq -70^{\circ}\text{C}$  for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature.

Avoid repeated freeze/thaw cycles.

## Contact



[www.abclonal.com](http://www.abclonal.com)