

RP10226K

Leader in Biomolecular Solutions for Life Science



# 26S Proteasome-mediated Protein Degradation Kit

Catalog No.: RP10226K

## Sequence Information

Species	Gene ID	Swiss Prot
Bovine		
proteins		

## Tags

## Synonyms

## Product Information

Source	Purification
Bovine red blood cells	Enzyme >95% by SDS-PAGE

## Endotoxin

## Formulation

Enzyme stocked in 20 mM Tris, 150 mM NaCl, 2 mM  $\beta$ ME, 10% Glycerol

## Reconstitution

## Background

This kit is designed for assaying 26S proteasome-mediated degradation of ubiquitinated proteins in vitro. It contains highly purified active bovine 26S proteasome that has been used by different research groups for assaying protein degradation.

## Basic Information

### Description

### Bio-Activity

### Storage

This product is stable at  $\leq -70^{\circ}\text{C}$  for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature.

Avoid repeated freeze/thaw cycles.

## Contact



[www.abclonal.com](http://www.abclonal.com)