

RP10238K

Leader in Biomolecular Solutions for Life Science



# Caspase-8, 10 Activity Fluorometric Assay Kit

Catalog No.: RP10238K

## Sequence Information

**Species**    **Gene ID**    **Swiss Prot**

**Tags**

**Synonyms**

## Background

Caspase-8, 10 Activity Fluorometric Assay Kit can be used for assaying caspase-8, 10, 3, 7 activities in cell/tissue extracts in a 96-well plate format. Each of the supplied fluorogenic substrate and the pan caspase inhibitor Z-VAD-FMK are sufficient for 250 x 100 µl reactions. The supplied TRAIL-treated cell lysate is used as a positive control, sufficient for 5 X 100 µl reactions.

## Basic Information

**Description**

**Bio-Activity**

**Storage**

This product is stable at  $\leq -70^{\circ}\text{C}$  for up to 6 months from the date of receipt. <br/>For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature.

Avoid repeated freeze/thaw cycles.

## Product Information

**Source**

Cell Lysates

**Purification**

Inhibitor and  
substrate > 98% by  
HPLC and NMR

**Endotoxin**

**Formulation**

**Reconstitution**

## Contact



[www.abclonal.com](http://www.abclonal.com)